

RAPIDCONTOUR



The Rapid Contour is a single-technology aesthetic device based on advanced Radio-Frequency (RF) technology. RF has been used in the aesthetic field for almost 20 years, and has proven itself to be the gold standard solution for a wide variety of conditions. SharpLight has fine-tuned the capabilities of RF and has developed a multi-polar head piece configuration that yields highly effective treatments that are completely pain-free.

Our Smart RF Hand piece studies your treatment impudence and adjusts the treatment time to be customized to your treatment method.

The Rapid Contour offers treatments are:

- Fast, pain-free and pleasant
- ▶ Effective, with visible results from the first treatment
- Lunch-time' procedures with no downtime
- ▶ 'Color blind' effective on all skin types
- Non-invasive

The Rapid Contour is equipped with two handpiecesthe Rapid Face and Rapid Body Our smart RF hand pieces study your treatment impudence and adjusts the treatment time according to your patient's current skin condition. The system has a state-of-the-art graphic user interface which shortens the learning curve and lets you simply set the treatment protocol . Designed to easily integrate into your practice, the Rapid Contour has a chic modern design and a small footprint.

Rapidly grow your practice

The business concept behind the Rapid Contour is to create a single technology device designated for wide body and face aesthetic treatments. As opposed to multi-technology platforms that require a greater financial investment, the single-technology alternative is more affordable. The Rapid Contour thus caters to practices that want lower entry price products or practices that want to add a stand-alone RF device to their existing aesthetic technologies. A key advantage of the Rapid Contour is that you will be utilizing all of its capabilities and no technology will go unused (which is not the case in multi-technology platforms).

With the Rapid Contour you can expect a fast return on investments. By offering a wide range of facial and body contouring solutions to your existing customers and also attract new customers to your practice. Treatments are fast, and no consumables are required.

About Radio-Frequency

Radio-frequency (RF) is a safe and effective technology used for a wide variety of aesthetic conditions. The advantage of RF is its ability to deliver targeted heat to deep tissue based on the tissue electrical resistance level properties.

Traditional bi-polar RF uses two polars through which the RF energy flows. SharpLight has designed an advanced multi-polar RF configuration in which the energy is equally distributed between. The advantage of the 6 poles configuration is that less but controlled energy is emitted. Resulting in an effective and comfortable treatment.

Skin rejuvenation:

RF can be used for a wide variety of indications such as skin rejuvenation, wrinkles, fine lines reduction, skin tightening & lifting and stretch marks. During treatment, the RF gently heats the dermis layer causing a contraction of the collagen fibers, which results in an immediate skin tightening effect. Concurrently, the heat stimulates the fibroblasts in the lower dermis to create neocollageogenesis and elastin, which, after repeated treatments, restores the skin's firmness and elasticity.

Circumferential reduction:

The body has a set number of fat cells, whereby weight fluctuation is a function of the fat cells' size. During treatment, RF energy heats up the sub-cutaneous adipose tissue. The heat increases the adipocyte metabolism which shrinks the fat tissue secreting fatty acids (liquid fat) which are then naturally disposed of via the body's blood and lymph systems. As the fat cells shrink, circumferential reduction is achieved.



"I own a small aesthetic practice and most of my business comes from body shaping. I chose the Rapid Contour since treatment with RF is pleasant and comfortable and doesn't require any consumables. My patients have been delighted with the results and my business turnover has increased dramatically."


Orly Zibli, Home of cosmetics and beauty, Israel



Cellulite is caused by the aggregation of enlarged adipocytes in the dermis layer and a weakening of the connective tissue. During RF treatment, energy heats the connective tissue that hold the cell structure and cause the cellulite and brakes them, simultaneously reducing the volume of the fat cell. Simultaneously, the RF energy strengthens the connective tissue. As a result, the skin's dimpling effect is treated, and it regains its smooth texture.

Cellulite (Human Skin)





Healthy skin

Skin with cellulite

Body contouring





Before

After



Number of treatments:11


Circumferential reduction







Cellulite





Number of treatments:7

SPECIFICATIONS

▶ Head options:

Face: 2 Flat electrodes

Body: 6 cylindrical electrodes 12mm 6 cylindrical electrodes 14mm

Power: up to 50 watt

ABOUT SHARPLIGHT

SharpLight Technologies Ltd. designs and manufactures non-invasive aesthetic medical solutions for treatment of a broad range of applications such as Hair Removal, Acne, Pigmented lesions, Vascular lesions, tattoo removal. The company is an expert in the aesthetic field and its proprietary DPC™ and VermaDerm™ technologies have redefined treatment effectiveness, comfort and speed. SharpLight offers a full line of stand-alone and multi-technology aesthetic workstations powered by DPC, RF, VermaDerm and Laser technology.

For more information visit SharpLightTech.com





SharpLight Technologies Ltd.
33 Lazarov St. Rishon Le Zion 7565435, Israel
T: +972.3.961.1969 | F: +972.3.952.8890 | E: info@sharplightech.com